





Benefits of implementing Powered Bed & Stretcher Movers in hospitals: the East Kent Hospitals story

Sharon Rindsland

Board of directors of the National Back Exchange & Trust Wide Moving and Handling Senior Co-Ordinator at East Kent Hospitals University Foundation NHS Trust. CEO of RQS Recognising Quality Solutions



AGENDA



- Hello and welcome
- Presentation from Sharon
- A quick overview of the Felgains Bed & Stretcher Mover range
- Live Q&A session

OBJECTIVES



- Musculoskeletal injuries from pushing and pulling
- The manual handling risks of manually pushing and pulling beds and stretchers.
- The impact this has on staff injury and staff turnover, and the subsequent financial and resource impacts.
- The benefits East Kent Hospitals have seen through implementing bed movers

HEALTH AND SAFETY STATISTICS



Key Figures for Great Britain (2021/22)

- 477,000 workers suffering from a work-related musculoskeletal disorder
- 61,713 injuries to employees reported under RIDDOR
- 36.8 million working days lost due to work-related illness and workplace injury
- £18.8 billion estimated cost of injuries and ill health from current working conditions (2019/20)



WHAT ARE MUSCULOSKELETAL DISORDERS (MSDS)?



Manual Handling:

• It is estimated that 1 in 3 accidents at work are caused by manual handling.

Acute:

 A sudden or one-time traumatic event or incident, e.g., slip, trip, fall or car wreck


Chronic or Cumulative:

- Injuries that occur over a <u>period of time</u> (months/years)
- & are caused by a combination of risk factors

MSDs affect ligaments, muscles, tendons, cartilage, blood vessels & nerves & spinal discs

PRIMARY RISK FACTORS FOR MSDS





THE CUMULATIVE EFFECT FROM PUSHING AND PULLING BEDS





Continued Exposure to Risk Factors



Time

EAST KENT HOSPITALS STORY (EKHUFT)



- 5 Hospital sites
- 1,400 bed spaces
- 10,500 staff
- All sites have slopes and tight corridors - Nightingale style




POOR ERGONOMIC DESIGN OF THE ENVIRONMENT OR THE EQUIPMENT BEING USED.



EKHUFT DEMANDS

- The job may be unnecessarily more difficult if it is not correctly and ergonomically matched to the employee;
- Staff would include manually pushing patient trolley or bed frequently over long distances or up and down slopes.





HOSPITAL BEDS



- The average bed/stretcher weighs 150kg
- Average bed weight of patients is 78kg

This means on average staff are pushing 228kg, this in on flat surface.

• As we have 1,400 beds, staff can be moving over a 100 beds a day.

Poor Posture






Potential Rotator cuff injury. Elbow away from the body

Arched back with added pulling forces lead to compression of the discs

Push and Pull Readings with Continuous movement




Force measured in NEWTONS



Start and Stop Pushing and Pulling Readings



Force measured in NEWTONS



RISK ASSESSMENT



- Staff should be involved in risk assessments
- Staff should also be included in the trial and selection of equipment that is purchased to reduce musculoskeletal risks.

MANUAL HANDLING GENERIC RISK ASSESSMENT

Task List <u>Number_EME</u> 1 Moving <u>Empty_bed</u> /trolley to go for repair	Is an assessment needed? to (is there a risk of <u>injury)</u> Yes
EME	

Task List Number EME 1 Moving Empty bed/trolley to go for repair to EME	Is an assessment needed? (is there a risk of <u>injury)</u> Yes
If 'Yes' continue. If 'No' the assessment need go	no further
Location: EME Trust Wide Personnel involved: EME Staff Date of Assessment: 08/11/13	Risk rating before changes 3L x 4C = 12 E Risk rating after changes 2L x 4C = 8 H If a powered mover used 1L x 2C= 2L Entered on Risk Register Yes
Operation covered by this assessment (detailed Description)	Diagrams, photographs, other information
straight sections of the route push or pull forces need to be applied at the foot end for tight turns. Free brakes are required for sideways manoeuvres There is high risk of injury if staff walk in	Table 2 shows the approximate increase in push forces that can be expected per 100 kg of load, on different slope angles. Table 2 Effect of slope angle on push force Slope gradient Push force (kg) increase (degrees) per 100 kg of laden trolley weight 1 2 3 5 9
front of moving beds/trolleys. Never walk in front with arms extended backwards to grip the foot board	7 12 As weight limit is above HSE guide, push beds/trolleys in straight lines and steer from
Many obstacles such as doorways, other equipment and human beings with restricted movement are present within the clinical settings and even more so in the public	the side using push/pull actions Do not walk in front of moving equipment or extend arms backwards
corridors	Ask for help if vision is obstructed or effort

Damage is caused to expensive equipment when the power plug is not removed or power

Excerpt from HSE Guidelines

cords wrapped around profiling sections of the equipment that then crush/cut through wiring

Another risk from pushing/pulling on a slope is

Ask for help if vision is obstructed or effort exceeds personal ability. Consider a powered mover for the task.

Prepare patient, attachments and power supply. Clear route of obstacles and clutter prior to task open doors Report, label and remove beds/trolleys that are

LOCa	List No: EME 1 Moving empty beartrolley between departments				
	Questions to be considered	Yes	No	Problems occurring from the task	Possible remedial action
	Insufficient rest or recovery?		1		
	Holding loads away from the trunk?	1		Risk of musculo-skeletal injury to lower back and	Manual Handling of Loads training for all staff plus
	Twisting?	V		shoulders due to movement and steering of a long load.	practical training for any powered movers
	Stooping?	1		Moving a poorly maintained bed/ trolley.	
	Reaching upwards?		٧	Chance to run lower limbs over if walking in front of	Report faulty beds/ trolleys to Estates/EME
	Large vertical Movement?		٧	moving equipment	
	Long carrying distance?		1	Minimum districts in a series and the series and the series and the series and the series are seri	Adjust height of equipment to enable clear vision
¥	Strenuous pulling or pushing?	1		Misaligned joints in spine, neck, shoulder and upper limbs if pulling the trolley	and naturally aligned joints
TASK	Unpredictable movement of loads?		1	innos ii puiling the trolley	Do not walk in front of moving bed/trolley. If
Ε.	Repetitive handling?		1	Obstacles to manoeuvre around especially humans	steering round tight corners push/pull from side
	Work rate imposed by process?	٧		Adverse cambers, door edging or clutter on floors	steering round tight corners pash/pail from side
	Jerking movement?	1		Adverse cambers, addressing or datter on noors	Use a powered mover if available. 2 staff required
	Static supporting for more than 2-3 seconds?		٧		if beyond an individuals physical capacity or
					personal height creates visual difficulties
					Clear area prior to task
					Set adequate time scales for task
	Loads are / do they				
	Heavy?	1			
ð	Unwieldy?	1		Risk of musculo-skeletal injury to lower back and	Powered mover or 2 staff required if beyond an
LOAD	Unstable / unpredictable?		٧	shoulders due to movement and steering of a long load.	individuals physical capacity or personal height
_	Intrinsically harmful?		1	Moving a poorly maintained bed/ trolley	creates visual difficulties Report faulty equipment to EME/Estates
	Require mechanical assistance?		V		Report faulty equipment to EME/Estates
	Are there				
	Constraints on posture?	V			Unplug power leads and secure cable to prevent
	Poor floors?	,	V	Possibility of twisting whilst moving long load in tight	entrapment in moving parts
눌	Variations in levels?		V	spaces	Clear the area of obstacles before undertaking
Ē	Hot / cold /. Humid conditions?		V		task. Open doors this may require a second
ź	Strong air movement?		v		person.
80	Poor lighting?		V	Obstacles/clutter or doorways	Keep shoulders and pelvis aligned when moving
ENVIRONMENT	Tripping hazard?	1			the trolley do not reach backwards
E				Human obstacles	Use pushing actions from rear or side to achieve
	Door Ho int				this
	Does the job		J		Una national matter
	Require unusual physical capability?	J	γ	May noce a rick to those who are prognant or have a	Use powered movers
	Hazardous to those with medical conditions / previous injury?	1	$\vdash \vdash$	May pose a risk to those who are pregnant or have a previous musculo-skeletal injury/health condition	Managere' programmy accomment required. Avoid
M	Hazardous to those who are pregnant?	γ	V	previous must my societal injurymeatin conducti	Managers' pregnancy assessment required. <u>Avoid</u> task if pregnant or ask for help
DO	Call for special information / training?		٧		task ii programi or ask for ficip
INDIVIDUAL					Work within your personal limits and report
ð					difficulties to Manager/Occupational Health
=					
Othe	Other factors: is movement or posture hindered by clothing or Personal Protective Equipment? No Could activity lead to cumulative strain: yes				

RISK RATINGS



LOW RISK

No immediate action. Update of the risk assessment to indicate the assessment outcome no further action

MEDIUM RISK Some level of assistance and/or equipment is required. State the assistance needed and prompt the assessor to identify the equipment (if required).

HIGH RISK

State the equipment required for many if not all activities.

BUSINESS CASE FOR EQUIPMENT



- Collect all information on cost
- Collect information on how many staff are needed to move patient
- Mitigate the risk, by weighing the risk

to staff and patient.



Level of Risk



Cost, time, trouble etc.


Occupational Health Referral



Date Typed: 27/06/2017 Clinic Date: 27/06/2017 Date Sent: 27/06/2017

Private & Confidential

WHIL



Dear

Re:

I saw for review in the Occupational Health Department at the William Harvey Hospital today. As you are aware was off work from September 2016 with back pain after moving a bed at work. MRI scanning showed him to have a prolapsed disc plus extensive degenerative change in his spine. He was seen by an orthopaedic surgeon at Medway Maritime Hospital recently who advised against heavy lifting tasks as a long term arrangement. The orthopaedic surgeon is making some changes to medication and is seeing him again for review on 15th August 2017. At that point they will consider offering medication treatment if indicated.

Many thanks for supporting return to work. He returned after I last met with him on 31st January 2017, and has been performing limited duties within the last met with him on 31st. In your email of 22nd February 2017 you informed me that the modifications.

would have performed a greater number of useful tasks more effectively. As I am sure you are aware under the Provision & Use of Work Equipment Regulations there is a legal responsibility on the

MANUAL HANDLING REPORT MOVEMENT OF BEDS THROUGHOUT WILLIAM HARVEY HOSPITAL



Ward/Department	William Harvey Hospital
Conducted by	Sharon Rindsland
	Moving and Handling Senior Co-Ordinator
Actioned by	Sharon Rindsland
	Moving and Handling Senior Co-Ordinator
Copies For	Incident reporting lead
	Nursing and portering leads
	Electric and Medical Engineering leads
	Legal Team
	Budget Directors
	Finances leads
Date 21/10/2016	29/10/2015

	CONTENT
Section 1	Introduction to Report
Section 2	General Comments
Section 3	Summary of Recommendations

ACCIDENTS AND LEGAL CASES



- Over a 12 month period we had 7 incident reports including one RIDDOR relating to staff injury pushing and pulling beds
- Risk Reduction staff must move the bed and patient with two staff+ depending on the size of patient and the bed.
- 3 Legal cases of staff taking action against the trust, minimum pay out 12,000 per person.
- 3 Impingement injuries within the shoulders

AVERAGE STAFF WAGES



- An average pay scale for a porter or a Health care assistant as it would be these staff moving the beds.
- Band 2: £21,406 (£10.95 an hour)
- To move the beds you need to release at least a minimum of 2 staff: £42,812



BALANCING THE COST OF EQUIPMENT WITH INJURY



- Time balancing cost
- Cost how many staff
 needed to undertake the
 task safely
- Benefits of Equipment up against Cost of injury



NO RISK ASSESSMENT = NO RISK

FELGAINS GZ10SL BED MOVER





- SWL: 500kg
- Slimline body allows each access into lifts and tight spaces.
- You don't need to change claws to attach to different beds.

FELGAINS GZS BED MOVER





- SWL: 600kg
- Fits around the bed frame, so no underbed clearance is needed.

FELGAINS GZ BARI BED MOVER





- SWL: 800kg
- For moving bariatric and heavier patients.

PURCHASE OF BED MOVERS



OBSTACLES

- Funding given bed movers purchased
- Porters wanted to use two staff or push Beds on their own.
- Competency training for bed movers completed



PURCHASE OF BED MOVERS - SUCCESS



- We now have 40 bed movers over 3 acute trusts
- Clinical staff and porters all now complete a incident report when their bed mover is not working



EKHUFT MANUAL HANDLING TEAM





Lee Sherwood



A quick overview of the Felgains Bed & Stretcher Movers



Live Q&A Session

REFERENCE LIST



- Eurostat. European Statistics on Accidents at Work (ESAW).
- Briggs, A.M.; Woolf, A.D.; Dreinhöfer, K.; Homb, N.; Hoy, D.G.; Kopansky-Giles, D.; March, L. Reducing the global burden of musculoskeletal conditions. Bull. World Health Organ. 2018, 96, 366.
- https://www.hse.gov.uk/statistics/causdis/msd.pdf
- https://www.hse.gov.uk/statistics/
- Serranheira, F.; Sousa-Uva, M.; Sousa-Uva, A. Hospital nurses' tasks and work-related musculoskeletal disorders symptoms: A detailed analysis. Work 2015, 51, 401–409
- Al-Eisawi K.W. Kerk C. Congleton J. Amendola A. Jenkins O. Gaines W. Factors Affecting Minimum Push and Pull Forces of Manual Carts Applied Ergonomics 30 (1999) P235-245
- Brace, T (2005) The Dynamics of Pushing and Pulling in the Workplace Assessing and Treating the problem. AAOHN Journal May 2005 Vol53 No 5 P224-229
- Cookson K. Rush A, (2011) The Guide to the Handling of People A systems Approach 6th ed Chapter 12 Page 194 People handling for Bariatrics

 A systems Approach Edited by Jacqui Smith, Published by Back Care Teddington Middlesex
- Ferreira, J. Boocock M. Gray M (2004) Review of the risks associated with pushing and pulling heavy loads Research Report 228 HMSO ISBN 0
 7176 2845 0 Great Britain Parliament (2012) Health and Social Care Act London HMSO
- https://www.clinicalservicesjournal.com/story/11285/pushing-and-pulling-keeping-staff-healthy







THANKS FOR JOINING!

Benefits of implementing powered Bed & Stretcher Movers in hospitals: the East Kent Hospitals story

