



# WHEELDRIVE 2.0

GO FURTHER.




mpulse





# WHEELDRIVE GO FURTHER



The WheelDrive 2.0 is unlike any other product on the market! Power assist wheels, designed to make driving easier.



The **WheelDrive 2.0** is a unique power add on wheel, with dual rim technology. A combination of power assistance and continuous push wheel drive.





The outer rim delivers extra power to every push. – offering speeds of up to **15km/h**. With less effort required per push, you'll enjoy longer journeys.



The inner rim provides consistent power, offering speeds of up to **6km/h**; allowing you the choice between propelling or consistent drive for those up hills.



Easy to fit on almost any wheelchair - without any additional parts or changes to your wheelchair.



Intuitive - It's just like regular self-propelling.



Go where you want with less effort.



# DUAL RIM TECHNOLOGY

You can naturally switch between both rims, so you'll always get the support required. Its unique design means the **WheelDrive 2.0** offers two drive rims to overcome the challenges of fatigue throughout the day.



**ASSIST Rim**

Traditional power assist rim which provides extra assistance when pushing. This means you can go further and faster per push - reducing strain on your body and leaving you with energy for your day. The outer rim offers speeds of up to **15km/h**.



**DRIVE Rim**

Unique to the WheelDrive, the inner rim produces a continuous drive option for when you're feeling tired of pushing. This allows you to reserve your energy for the day ahead. The inner rim offers speeds of up to **6km/h**.





# INTUITIVE TO USE

## FUNCTIONAL BUTTONS

With an easily accessible power button, you can simply switch the system on and off and set your driving profile with three different speeds.

## ROLL BACK STOP

When propelling up a hill with the outer rim, the WheelDrive 2.0 enables a roll back stop. The roll back stop allows you to position your hands for the next push without rolling back; which will improve comfort and safety when traveling uphill.

## DOWNHILL SPEED LIMITER

The WheelDrive 2.0 provides a top speed cut out. When travelling down a hill, the WheelDrive keeps a consistent speed which keeps you in control – with the inner rim limited to **6 km/h** and the outer rim limited to **15 km/h**. This feature maintains your safety and confidence in downhill driving.



# EMPOWERED DESIGN

The control, functionality and looks of a manual wheelchair combined with leading power assist and drive technology. With the WheelDrive 2.0, you get the best of both worlds. Go further – and enjoy an effortless ride!







### // EASY TRANSPORT & AIR TRAVEL

The WheelDrive 2.0 can be quickly removed from your wheelchair for easy transportation. The WheelDrive also fulfills IATA requirements and may be carried as checked baggage. (Remove batteries before transport.)



### // SAFETY DURING TRANSPORT

This product has been successfully crash tested and fulfills the performance requirements for ISO 7176-19.



### // FITS ON ALMOST ANY WHEELCHAIR

The WheelDrive 2.0 can be easily mounted on virtually any lightweight or standard wheelchair, without any additional parts or changes to the chair.



### // SMART BATTERY TECHNOLOGY

The WheelDrive 2.0 is powered by two easily detachable lightweight batteries which deliver exceptional performance. These batteries can be fully charged in two hours. WheelDrive's detachable lightweight batteries allow the user to propel the WheelDrive manually, while the batteries are charging.

# TECHNICAL DATA



## WheelDrive // WHEELCHAIR POWER ASSIST

<b>Max. load capacity</b>	130 kg	<b>Max. speed forward</b>	10 km/h (assist rim); 6 km/h (drive rim)
<b>Total extra width compared to manual wheels (set)</b>	43 mm	<b>Max. range (ISO 7176-4)</b>	20 km (assist rim); 12 km (drive rim) (dependant on user weight, driving programme, weather and road conditions, condition of battery and tyres, etc.)
<b>Weight of batteries (set)</b>	3.6 kg		
<b>Total weight including battery (one wheel)</b>	11.5 kg	<b>Max. kerb climbing, driving backwards</b>	50 mm
<b>Transport weight excluding battery (one wheel)</b>	9.7 kg	<b>Charging time (full charge)</b>	Approx. 2 hours
<b>Wheel diameter</b>	24"	<b>Intended use and environment</b>	For people with limited mobility; self-propelling; indoor and outdoor use. Wheelchair type class B
<b>Crash tested</b>	Yes	<b>Reimbursement Code</b>	18.99.05.2001



### Sunrise Medical Limited

Thorns Road, Brierley Hill, West Midlands, DY5 2LD  
 Tel.: +44 (0) 845 605 66 88  
 Fax: +44 (0) 845 605 66 89  
 E-Mail: [enquires@sunmed.co.uk](mailto:enquires@sunmed.co.uk)  
[www.sunrisemedical.co.uk](http://www.sunrisemedical.co.uk)



For further information on the full specification, options and accessories, please refer to the order form.  
 All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.co.uk](http://www.SunriseMedical.co.uk)