


# Pallas slings



### Universal sling

Support level: Medium (supportive of trunk and thighs)

Positioning: Upright or standing

Use: Standing or seated to seated transfers

Characteristics of resident:

Some weight bearing ability with low trunk control

Attachment: Loop

SWL: 200kg/31 stone



#### Chest sling

Support level: Low (supportive of trunk)

Positioning: Upright or standing


Use: Standing or seated to seated transfers

Characteristics of resident:

Some weight bearing ability with good trunk control

Attachment: Loop

SWL: 200kg/31 stone



### Fleece sling

Used in the same way as the Chest sling, but provides additional comfort for the resident.

Support level: Low (supportive of trunk)

Positioning: Upright or standing

Use: Standing or seated to seated transfers

Characteristics of resident:

Some weight bearing ability with good trunk control

Attachment: Loop

SWL: 200kg/31 stone

## Pallas sling models

Universal sling		
FEL30-234	Medium	Polyester fabric
FEL30-113	Large	Polyester fabric
Chest sling		
FEL30-074	Small	Spacer fabric
FEL30-075	Medium	Spacer fabric
FEL30-076	Large	Spacer fabric
FEL30-077	Extra Large	Spacer fabric
Fleece sling		
FEL30-070	Small	Fleece fabric
FEL30-071	Medium	Fleece fabric
FEL30-072	Large	Fleece fabric
FEL30-073	Extra Large	Fleece fabric

#### Colour coding

Small	
Medium	
Large	
Extra Large	

## Have any questions?

Call us on 01473 741144 or send an email to customerservices@felgains.com.